



UZ
LEUVEN

HK VAATHEELKUNDE



Leuven, 23 September 2011

Dear,

EVRF is a wonderful technique for treating veins insufficiency. The Saphena Magna and Parva, as well as smaller reticular veins. This is a great advantage, since the correct amount of energy can be given per diameter to obtain the occlusion of the vein. And also shorter segments can be treated, as the tip of the catheter is only 0,5cm long. The catheter (CR45i) for the treatment of stamvarices is very smooth and easy to carry, even over a twist. The catheter lies good in de hand and with the roller system it's easy to control. The dozens of patient treated here, were all satisfied with the result. We will soon launch a pilot study with 40 VSM's. The focus will be on the Quality of Life score, pain score and occlusion rate.

Yours sincerely,

Dr. Sarah Thomis

Flebologe, Vascular Surgery Service

UZ Leuven